

The <b>GOODELL</b> laboratory		
<b>Author</b>	Chris Benton	May 8, 2009
<b>Title</b>	<b>Bone marrow isolation, crushing technique</b>	
<b>Introduction</b>	This protocol describes the isolation of whole bone marrow from whole mice using a bone crushing technique. This protocol may be particularly useful to isolate large numbers of bone marrow cells, or to isolate bone marrow from mice for which a bone flushing technique is difficult [i.e. young or small mice].	
<b>Materials</b>	<ol style="list-style-type: none"> <li>1. Dissecting tools, including scissors, various sized tweezers.</li> <li>2. Mortar and pestle, sterilized.</li> <li>3. HBSS+: Hanks Balanced Salt Solution (from Gibco) + 2% Heat-inactivated Fetal Bovine Serum + 10mM Hepes + Pen/Strep.</li> <li>4. PBS and PBS+: PBS + 2% Heat-inactivated Fetal Bovine Serum.</li> <li>5. Collagenase/Dispase [Roche cat# 11097113001].</li> <li>6. 50mL conical centrifuge tubes.</li> <li>7. 40 micron sterile cell strainers.</li> <li>8. Table top centrifuge capable of spinning falcon tubes at ~500xg.</li> <li>9. Tissue culture dishes.</li> </ol>	
<b>Protocol</b>		<i>Notes</i>
1.	Prepare bone marrow from mice starting by dissecting tibias, femurs and iliac crests, cleaned of all muscle and connective tissue. Place intact bones into a tissue culture dish on ice with PBS+. Continue by dissecting the arms and scapulas, again cleaning away all muscle and connective tissue. Dissect sternum and spine away from ribs, ensuring spine is also devoid of spinal cord.	<i>We do not typically dissect ribs, however any bones maybe used in this protocol, so long as they are dissected free from all soft tissue.</i>
2.	Transfer cleaned bones to sterile mortar, containing 1mL PBS+ per 3 bones.	<i>After trisecting the spine, we find most mice give rise to ~18 bones, so we generally use 6mL PBS+ per mouse.</i>
3.	Using the pestle, crush bones into bone fragments,	



