

	<b>THE GOODELL LABORATORY</b>	
<b>Author</b>	Karen Lin	July 30, 2005
<b>Title</b>	<b>Production of MSCV Retrovirus</b>	

	<p><b>Dissolve DNA in OPTI medium, incubate for 5 minutes</b></p> <p>MSCV vector            2ug /RXN  Eco-pCL vector        2ug /RXN  OPTI medium            250ul/ RXN</p> <p>3. <u>Make LipofectAmine-OPTI mixture (Solution B)</u>  Dissolve LipofectAmine in OPTI medium  10ul LipofectAmine + 250 ul OPTI /RXN  <b>Incubate for 5 minutes.</b></p> <p>4. Mix DNA (Solution A) and LipofectAmine mixture (Solution B), incubate in RT for 20 minutes.</p> <p>5. Gently apply DNA-LipofectAmine mixture (500ul/RXN) onto the 293T cells in the 6 well plates. Avoid aspiration.</p>	
<p>3.</p>	<p><u>d 1 Replace medium</u>  Replace medium-containing medium, 2ml/ RXN</p>	<p>293T cells may be easily detached from the plate. Removed the old medium. Do not let cells sit without having new medium added so you will have to be quick. When adding new medium, gently tilt the plate, and add the medium to the wall of well so the medium drop would not disturb the cells.</p>
	<p>4. <u>d 2 Harvest virus</u>  48 hours after transfection, harvest supernatant. To exclude cell debris, one can 1) spin supernatant in 4°C, 2000rpm for 10 minutes, and transfer the supernatant to the freezing tubes; or 2) filter the virus through a 0.45uM syringe filter.</p>	<p>To prevent losing viral titer in the future during viral titering and infection experiment, repeated thaw-freeze cycles need to be</p>